


Acumar™

Range of Motion

Measurement Illustrations



ACUMAR
By



Lafayette Instrument Company
3700 Sagamore Parkway North
Lafayette, IN 47904
(765) 423-1505 FAX: (765) 423-4111
www.lafayetteinstrument.com

Measurement Sections

The Spine:

Cervical:	Page
	3
Flexion, Extension	
Left, Right Lateral	
Left, Right Rotation	
Thoracic:	4
Flexion, Extension	
Left, Right Lateral	
Left, Right Rotation	
Lumbar:	5
Flexion, Extension	
Left, Right Lateral	

Upper Extremities:

Shoulder:	6-7
Flexion, Extension	
Abduction, Adduction	
External, Internal Rotation	
Elbow:	8
Flexion, Extension	
Wrist:	9-10
Flexion, Extension	
Radial, Ulnar Deviation	

Lower Extremities

Hip:	11-12
Flexion, Extension	
Abduction, Adduction	
External, Internal Rotation	
Knee:	13
Flexion, Extension	
Ankle:	14-15
Flexion, Extension	
Inversion, Eversion	



Disclaimer





This manual is designed solely to illustrate the use of the Acumar Digital Inclinometer, and is not intended nor implied to describe medical examination procedures. Selecting methods of examination and interpreting measured results are the responsibility of practitioners.

Acumar will not warrant that any procedures described herein meet any current or past medical examination guides or procedures required by state or any local governments or other organizations. Acumar expressly disclaims all liability arising from use of this manual. Acumar Digital Inclinometer or related product(s) warranties are as stated in the warranty in the latter part of the users manual.

The use of Acumar Digital Inclinometers to take objective measurements combined with the AccuCapture™ excel template can produce patient files. Like any other patient data, if protection of the privacy of patient files is required, it is the sole responsibility of the provider and of the practitioner. Acumar disclaims all liability for any and all claims.

Acumar™ does not express or imply that the reference values on these forms are accurate, relevant, or clinically significant and hereby disclaims all liability for their clinical use herein. Medical judgment should be the ultimate faculty of the practitioner.


The Spine: Cervical Range of Motion

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on the top of the head and then the companion unit on the T1 spinous process. Press the ZERO button and then ask the individual to FLEX forward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times.</p>		<p>Record data, then proceed to the next test. To start the next sequence of tests, press START (off), press START (on) to clear memory. With individual in neutral position, press ZERO. Then ask the individual to EXTEND backward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times</p>
		
Left Lateral	Neutral	Right Lateral
<p>While individual is in the neutral position place the main unit on the top of the head and then the companion unit (hidden) on the T1 spinous process. Press the ZERO button and then ask the individual to FLEX to the LEFT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the top of the head and then the companion unit (hidden) on the T1 spinous process. Press the ZERO button and then ask the individual to FLEX to the RIGHT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>
		
Left Rotation	Neutral Position	Right Rotation
<p>While individual is in the neutral supine position on a flat exam table, place the main unit on the forehead. Press the ZERO button and then ask the individual to ROTATE his or her head to the LEFT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. <i>Measuring cervical rotation requires stabilization of the trunk in the supine position with the shoulders on the table.</i></p>		<p>While individual is in the neutral supine position on a flat exam table, place the main unit on the forehead. Press the ZERO button and then ask the individual to ROTATE his or her head to the RIGHT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. <i>Measuring cervical rotation requires stabilization of the trunk in the supine position with the shoulders on the table.</i></p>

The Spine: Thoracic Range of Motion

					
Flexion		Neutral		Extension	
<p>While individual is in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to FLEX forward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times.</p>			<p>While individual is in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to EXTEND backward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times.</p>		
					
Left Lateral				Right Lateral	
<p>While individual is in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to FLEX to the LEFT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to FLEX to the RIGHT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>		
					
Left Rotation				Right Rotation	
<p>While individual is seated or standing and in a forward flexed position with the thoracic spine horizontal, in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to ROTATE to the LEFT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>			<p>While individual is seated or standing and in a forward flexed position with the thoracic spine horizontal, in the neutral position place the main unit on the T1 spinous process and then the companion unit on the T12 spinous process. Press the ZERO button and then ask the individual to ROTATE to the RIGHT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>		

The Spine: Lumbar Range of Motion

					
Flexion		Neutral		Extension	
<p>While individual is in the neutral position place the main unit on the T12 spinous process and then the companion unit on the S1 spinous process. Press the ZERO button and then ask the individual to FLEX forward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times.</p>			<p>While individual is in the neutral position place the main unit on the T12 spinous process and then the companion unit on the S1 spinous process. Press the ZERO button and then ask the individual to EXTEND backward completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data. Repeat two to five times.</p>		
					
Left Lateral		Neutral		Right Lateral	
<p>While individual is in the neutral position place the main unit on the T12 spinous process and then the companion unit on the S1 spinous process. Press the ZERO button and then ask the individual to FLEX to the LEFT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the T12 spinous process and then the companion unit on the S1 spinous process. Press the ZERO button and then ask the individual to FLEX to the RIGHT completely. When the individual has reached the maximum range of motion, press the HOLD button to store the data.</p>		

Upper Extremities: Shoulder Range of Motion: LEFT SHOULDER

Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to FLEX their LEFT SHOULDER completely. Press the HOLD button to store the data. Repeat two to five times</p>		<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to EXTEND their LEFT SHOULDER completely. Press the HOLD button to store the data. Repeat two to five times.</p>
Adduction	Neutral	Abduction
<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ADDUCT their LEFT SHOULDER completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ABDUCT their LEFT SHOULDER completely. Press the HOLD button to store the data.</p>
External Rotation	Neutral	Internal Rotation
<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ROTATE their LEFT SHOULDER EXTERNALLY completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the LEFT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ROTATE their LEFT SHOULDER INTERNALLY completely. Press the HOLD button to store the data.</p>

Upper Extremities: Shoulder Range of Motion: RIGHT SHOULDER

					
Flexion		Neutral		Extension	
<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to FLEX their RIGHT SHOULDER completely. Press the HOLD button to store the data. Repeat two to five times</p>			<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to EXTEND their RIGHT SHOULDER completely. Press the HOLD button to store the data. Repeat two to five times.</p>		
					
Adduction		Neutral		Abduction	
<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ADDUCT their RIGHT SHOULDER completely. Press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ABDUCT their RIGHT SHOULDER completely. Press the HOLD button to store the data.</p>		
					
External Rotation		Neutral		Internal Rotation	
<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ROTATE their RIGHT SHOULDER EXTERNALLY completely. Press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the RIGHT humerus. For single inclinometry stabilize movement across the shoulders, for dual inclinometry, place the companion unit on the clavicle. Press the ZERO button and then ask the individual to ROTATE their RIGHT SHOULDER INTERNALLY completely. Press the HOLD button to store the data.</p>		

Upper Extremities: Elbow Range of Motion: LEFT ELBOW

			
Flexion		Neutral	Extension
<p>While individual is in the neutral position place the main unit on the LEFT ulna or radius For single inclinometry stabilize movement in the humerus, for dual inclinometry, place the companion unit on the humerus Press the ZERO button and then ask the individual to FLEX the LEFT ELBOW completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on the LEFT ulna or radius For single inclinometry stabilize movement in the humerus, for dual inclinometry, place the companion unit on the humerus Press the ZERO button and then ask the individual to EXTEND the LEFT ELBOW completely. Press the HOLD button to store the data. Repeat two to five times.</p>	

Upper Extremities: Elbow Range of Motion: RIGHT ELBOW

			
Flexion		Neutral	Extension
<p>While individual is in the neutral position place the main unit on the RIGHT ulna or radius For single inclinometry stabilize movement in the humerus, for dual inclinometry, place the companion unit on the humerus Press the ZERO button and then ask the individual to FLEX the RIGHT ELBOW completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the RIGHT ulna or radius For single inclinometry stabilize movement in the humerus, for dual inclinometry, place the companion unit on the humerus Press the ZERO button and then ask the individual to EXTEND the RIGHT ELBOW completely. Press the HOLD button to store the data.</p>	

Upper Extremities: Wrist Range of Motion: LEFT WRIST

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on their LEFT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to FLEX the LEFT WRIST completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on their LEFT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to EXTEND the LEFT WRIST completely. Press the HOLD button to store the data. Repeat two to five times.</p>
		
Radial Deviation	Neutral	Ulnar Deviation
<p>While individual is in the neutral position place the main unit on their LEFT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to DEVIATE the LEFT RADIUS completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on their LEFT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to DEVIATE the LEFT ULNA completely. Press the HOLD button to store the data.</p>

Upper Extremities: Wrist Range of Motion: RIGHT WRIST

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on their RIGHT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to FLEX the RIGHT WRIST completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on their RIGHT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to EXTEND the RIGHT WRIST completely. Press the HOLD button to store the data. Repeat two to five times.</p>
		
Radial Deviation	Neutral	Ulnar Deviation
<p>While individual is in the neutral position place the main unit on their RIGHT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to DEVIATE the RIGHT RADIUS completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on their RIGHT phalanges. For single inclinometry stabilize movement in the forearm, for dual inclinometry, place the companion unit on the ulna or radius Press the ZERO button and then ask the individual to DEVIATE the RIGHT ULNA completely. Press the HOLD button to store the data.</p>

Lower Extremities: Hip: LEFT HIP

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to FLEX their LEFT HIP completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to EXTEND their LEFT HIP completely. Press the HOLD button to store the data. Repeat two to five times.</p>
		
Abduction	Neutral	Adduction
<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to ABDUCT their LEFT HIP completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to ADDUCT their LEFT HIP completely. Press the HOLD button to store the data.</p>
		
External Rotation	Neutral	Internal Rotation
<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to EXTERNALLY ROTATE their LEFT HIP completely. Press the HOLD button to store the data. Repeat two to five times. When you are ready to send data to the computer, press the send button while holding the main unit backside to the IR receiver.</p>		<p>While individual is in the neutral position place the main unit on the LEFT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to INTERNALLY ROTATE their LEFT HIP completely. Press the HOLD button to store the data. Repeat two to five times. When you are ready to send data to the computer, press the send button while holding the main unit backside to the IR receiver.</p>

Lower Extremities: Hip: RIGHT HIP

Flexion		Neutral	Extension	
<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to FLEX their RIGHT HIP completely. Press the HOLD button to store the data. Repeat two to five times.</p>			<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to EXTEND their RIGHT HIP completely. Press the HOLD button to store the data. Repeat two to five times.</p>	
Abduction		Neutral	Adduction	
<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to ABDUCT their RIGHT HIP completely. Press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to ADDUCT their RIGHT HIP completely. Press the HOLD button to store the data.</p>	
External Rotation		Neutral	Internal Rotation	
<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to EXTERNALLY ROTATE their RIGHT HIP completely. Press the HOLD button to store the data.</p>			<p>While individual is in the neutral position place the main unit on the RIGHT femur. Stabilize movement across the pelvis. Press the ZERO button and then ask the individual to INTERNALLY ROTATE their RIGHT HIP completely. Press the HOLD button to store the data.</p>	

Lower Extremities: Hip: LEFT KNEE

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on their LEFT TIBIA. For single inclinometry stabilize movement in the upper leg, for dual inclinometry, place the companion unit on the femur. Press the ZERO button and then ask the individual to FLEX the LEFT KNEE completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on their LEFT TIBIA. For single inclinometry stabilize movement in the upper leg, for dual inclinometry, place the companion unit on the femur. Press the ZERO button and then ask the individual to EXTEND the LEFT KNEE completely. Press the HOLD button to store the data. Repeat two to five times.</p>

Lower Extremities: Hip: RIGHT KNEE

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on their RIGHT TIBIA. For single inclinometry stabilize movement in the upper leg, for dual inclinometry, place the companion unit on the femur. Press the ZERO button and then ask the individual to FLEX the RIGHT KNEE completely. Press the HOLD button to store the data. Repeat two to five times. When you are ready to send data to the computer, press the send button while holding the main unit backside to the IR receiver.</p>		<p>While individual is in the neutral position place the main unit on their RIGHT TIBIA. For single inclinometry stabilize movement in the upper leg, for dual inclinometry, place the companion unit on the femur. Press the ZERO button and then ask the individual to EXTEND the RIGHT KNEE completely. Press the HOLD button to store the data. Repeat two to five times. When you are ready to send data to the computer, press the send button while holding the main unit backside to the IR receiver.</p>

Lower Extremities: Ankle: LEFT ANKLE

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on the bottom of their LEFT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to FLEX the LEFT ANKLE completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on the bottom of their LEFT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to EXTEND the LEFT ANKLE completely. Press the HOLD button to store the data.</p>
		
Inversion	Neutral	Eversion
<p>While individual is in the neutral position place the main unit on the bottom of their LEFT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to INVERT the LEFT ANKLE completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the bottom of their LEFT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to EVERT the LEFT ANKLE completely. Press the HOLD button to store the data.</p>

Lower Extremities: Ankle: RIGHT ANKLE

		
Flexion	Neutral	Extension
<p>While individual is in the neutral position place the main unit on the bottom of their RIGHT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to FLEX the RIGHT ANKLE completely. Press the HOLD button to store the data. Repeat two to five times.</p>		<p>While individual is in the neutral position place the main unit on the bottom of their RIGHT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to EXTEND the RIGHT ANKLE completely. Press the HOLD button to store the data. Repeat two to five times.</p>
		
Inversion	Neutral	Eversion
<p>While individual is in the neutral position place the main unit on the bottom of their RIGHT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to INVERT the RIGHT ANKLE completely. Press the HOLD button to store the data.</p>		<p>While individual is in the neutral position place the main unit on the bottom of their RIGHT FOOT. For single inclinometry stabilize movement in the lower leg, for dual inclinometry, place the companion unit on the tibia Press the ZERO button and then ask the individual to EVERT the RIGHT ANKLE completely. Press the HOLD button to store the data.</p>