Lafayette Pegboard User's Manual





P.O. Box 5729 Lafayette, IN 47903 USA

Tel: 765.423.1505 • 800.428.7545

Fax: 765.423.4111

info@lafayetteinstrument.com www.lafayetteinstrument.com

Lafayette Pegboard Test User's Manual

Description

The Lafayette Pegboard Test can be used as a manipulative dexterity test. It is designed to measure a timed unilateral motor response. This unit consists of 25 round holes on a board, in which subject is to place pegs as quickly as possible. This test requires visual-motor coordination, and speed should improve over successive trials.

Instructions for the Test

The pegboard is placed in mid-line with the subject so that the board is at the edge of the table and peg tray immediately above the board. The examiner explains the test:

"THIS IS A PEGBOARD AND THESE ARE THE PEGS." (The examiner points out each and then picks up one o the pegs and continues.)

"ALL THE PEGS ARETHE SAME. WHAT YOU MUST DO IS PUTTHESE PEGS INTO THE HOLES LIKETHIS." (The examiner demonstrates by filling the top row.)

(Remove the pegs, putting them back into the tray.)

"WHEN I SAY GO, BEGIN HERE AND PUT THE PEGS INTO THE BOARDS AS FAST AS YOU CAN, USING ONLY YOUR (DOMINANT) HAND. FILL THE TOP ROW COMPLETELY FROM THIS SIDE TO THIS SIDE. DO NOT SKIP ANY; FILL EACH ROW THE SAME WAY YOU FILLED THE TOP ROW. ANY QUESTIONS? READY, AS FAST AS YOU CAN. GO."

Supplementary Instructions

For the right hand trial, the examiner demonstrates that the pegs are placed from subject's left to right, and from right to left for the left-hand trial. The dominant hand trial is administered first, followed by the non-dominant hand trial.

The examiner encourages the subject to perform the task as quickly as possible. The pegs must be put in the board in exact order. Frequently, it will be necessary to point out the first hole of a new row, particularly during the non-dominant hand trial. Only one peg is to be picked up a time and the subject should immediately be told if more than one is picked up.

Also, only one hand is to be used. Occasionally, a subject will attempt to use his or her other hand to help. It may be necessary to tell the subject to keep the hand on his or her lap, or for the examiner to hold it. If necessary, the board should be held steady for the patient. In the case of sever motor impairment, the subject should attempt the task just to see if any of the pegs can be put in. Any factor which mat effect the subjects performance should be noted, e.g. sore finger, bandage, etc.

If a peg is dropped to the floor, the examiner should not make an attempt to pick it up during the trial; rather, one of the pegs correctly placed should be taken out and used again.

Model 32027 User's Manual

Scoring and Interpretation

Record, in seconds, the length of time required to fill an holes as the subject's score.

Test scores can be used in an individual evaluative or therapeutic selection, the ideal procedure is to establish its validity locally, by testing all newly hired employees and correlating scores with their subsequent performance (supervisor ratings or duration of employment). This approach to validation requires that test scores should not be used to select employees until evidence has accumulated of its validity, and that test scores be inaccessible to supervisors or others who affect the worker's ratings or job longevity. An alternative validation procedure is to administer the test to all present employees and correlate scores with ratings or with subsequent performance. (Anastasi, 1982, pp. 54-101). The Lafayette pegboard test should correlate most highly with those jobs, which require speed, finger dexterity and manual dexterity. It should be of relevance to performance on assembly and machine operating jobs.

Reference

Anastasi, Anne. Psychological Testing. (5th ed.). New York: Macmillan, 1982.



Terms and Conditions

LIC Worldwide Headquarters

Toll-Free: (800) 428-7545 (USA only) Phone: (765) 423-1505

Phone: (765) 423-150 Fax: (765) 423-4111

Email: sales@lafayetteinstrument.com export@lafayetteinstrument.com (Outside the USA)

Mailing Address:

Lafayette Instrument Company

PO Box 5729 Lafayette, IN 47903, USA

Lafayette Instrument Europe:

Phone: +44 1509 817700 Fax: +44 1509 817701

Email: eusales@lafayetteinstrument.com

Phone, Fax, Email or Mail-in Orders

All orders need to be accompanied by a hard copy of your purchase order. All orders must include the following information:

- Quantity
- Part Number
- Description
- Your purchase order number or method of pre-payment
- Your tax status (include tax-exempt numbers)
- Shipping address for this order
- · Billing address for the invoice we'll mail when this order is shipped
- Signature and typed name of person authorized to order these products
- Your telephone number
- Your email address
- Your FAX number

Domestic Terms

There is a \$50 minimum order. Open accounts can be extended to most recognized businesses. Net amount due 30 days from the date of shipment unless otherwise specified by us. Enclose payment with the order; charge with VISA, MasterCard, American Express, or pay COD. We must have a hard copy of your purchase order by mail, E-mail or fax. Students, individuals and private companies may call for a credit application.

International Payment Information

There is a \$50 minimum order. Payment must be made in advance by: draft drawn on a major US bank; wire transfers to our account; charge with VISA, MasterCard, American Express, or confirmed irrevocable letter of credit. Proforma invoices will be provided upon request.

Exports

If fordering instrumentation for use outside the USA, please specify the country of ultimate destination, as well as the power requirements (110V/60Hz or 220V/50Hz). Some model numbers for 220V/50Hz will have a "*C" suffix

Quotations

Quotations are supplied upon request. Written quotations will include the price of goods, cost of shipping and handling, if requested, and estimated delivery time frame. Quotations are good for 30 days, unless otherwise noted. Following that time, prices are subject to change and will be re-quoted at your request.

Cancellations

Orders for custom products, custom assemblies or instruments built to customer specifications will be subject to a cancellation penalty of 100%. Payment for up to 100% of the invoice value of custom products may be required in advance. Cancellation for a standard Lafayette Instrument manufactured product once the product has been shipped will normally be assessed a charge of 25% of the invoice value, plus shipping charges. Resell items, like custom products, will be subject to a cancellation penalty of 100%.

Exchanges and Refunds

Please see the cancellation penalty as described above. No item may be returned without prior authorization of Lafayette Instrument Company and a completed Return Form. A copy of the Return Form or your assigned Return # (you will receive this via email after submitting the form) must be included with the returned goods. The merchandise should be packed well and fully insured. Unopened merchandise may be returned prepaid within thirty (30) days after receipt of the item and in the original shipping carton. Collect shipments will not be accepted. Returned products must be in saleable condition, and credit is subject to inspection of the merchandise.

Repair

Instrumentation may not be returned without prior authorization by Lafayette Instrument Company and a completed Return Form. When you

complete the Form, or call Lafayette Instrument, you will receive a Return #. Your Return # number will be good for 30 days. Address the shipment to:

Lafayette Instrument Company 3700 Sagamore Parkway North Lafayette, IN 47904, USA.

Shipments cannot be received at the LIC PO Box. Items should be packed well, insured for full value, and returned along with a copy of the Return Form or the Return #. An estimate of repair will be given prior to completion ONLY if requested in an enclosed cover letter. We must have a completed purchase order by mail or fax, or repair work cannot commence for non-warranty repairs.

Damaged Goods

Damaged instrumentation should not be returned to Lafayette Instrument prior to a thorough inspection. If a shipment arrives damaged, note damage or delivery bill and have the driver sign it to acknowledge the damage. Contact the delivery service, and they will file an insurance claim. If damage is not detected at the time of delivery, contact the carrier/shipper and request an inspection within 10 days of the original delivery. Please call the Lafayette Instrument Customer Service Department for replair or replacement of the damaged merchandise.

Limited Warranty

Lafayette Instrument Company warrants equipment manufactured by the company to be free of defects in material and workmanship for a period of one year from the date of shipment, except as provided hereinafter. The original manufacturer's warranty will be honored by Lafayette Instrument for items not manufactured by Lafayette Instrument Company, i.e. resell items. This assumes normal usage under commonly accepted operating parameters and excludes consumable products.

Warranty period for repairs or used instrumentation purchased from Lafayette Instrument is 90 days. Lafayette Instrument Company agrees either to repair or replace, at its sole option and free of part charges to the customer, instrumentation which, under proper and normal conditions of use, proves to be defective within the warranty period. Warranty for any parts of such repaired or replaced instrumentation shall be covered under the same limited warranty and shall have a warranty period of 90 days from the date of shipment or the remainder of the original warranty period whichever is greater. This warranty and remedy are given expressly and in lieu of all other warranties, expressed or implied, of merchantability or fitness for a particular purpose and constitutes the only warranty made by Lafayette Instrument Company.

Lafayette Instrument Company neither assumes nor authorizes any person to assume for it any other liability in connection with the sale, installation, service or use of its instrumentation. Lafayette Instrument Company shall have no liability whatsoever for special, consequential, or punitive damages of any kind from any cause arising out of the sale, installation, service or use of its instrumentation. All products manufactured by Lafayette Instrument Company are tested and inspected prior to shipment. Upon prompt notification by the Customer, Lafayette Instrument Company will correct any defect in warranted equipment of its manufacture either, at its option, by return of the item to the factory, or shipment of a repaired or replacement part. Lafayette Instrument Company will not be obliged, however, to replace or repair any piece of equipment, which has been abused, improperly installed, altered, damaged, or repaired by others. Defects in equipment do not include decomposition, wear, or damage by chemical action or corrosion, or damage incurred during shipment.

Limited Obligations Covered by this Warranty

- In the case of instruments not of Lafayette Instrument Company manufacture, the original manufacturer's warranty applies.
- Shipping charges under warranty are covered only in one direction. The customer is responsible for shipping charges to the factory if return of the part is required.
- This warranty does not cover damage to components due to improper installation by the customer.
- Consumable and or expendable items, including but not limited to electrodes, lights, batteries, fuses, O-rings, gaskets, and tubing, are excluded from warranty.
- Failure by the customer to perform normal and reasonable maintenance on instruments will void warranty claims.
- If the original invoice for the instrument is issued to a company that is not the company of the end user, and not an authorized Lafayette Instrument Company distributor, then all requests for warranty must be processed through the company that sold the product to the end user, and not directly to Lafayette Instrument Company.

Export License

The U.S. Department of Commerce requires an export license for any polygraph system shipment with an ULTIMATE destination other than: Australia, Japan, New Zealand or any NATO Member Countries. It is against U.S. law to ship a Polygraph system to any other country without an export license. If the ultimate destination is not one of the above listed countries, contact us for the required license application forms.